


# NX NAIL TECHNICAL SPECIFICATIONS

**FIELD** ORTHOPAEDICS



Nail Size	Indication	Proximal Head $\varnothing$	Shaft $\varnothing$	Distal Tip $\varnothing$	Shortest Implant	Longest Implant	K-Wire $\varnothing$	Drill $\varnothing$	Screw Driver
2.0	Phalanx, Scaphoid	2.7	1.5	2.0	12	40	0.8	1.5	1.5mm Hex
2.5	Phalanx, Scaphoid, Metacarpal	3.15	1.75	2.5	12	60	1.0	1.75	1.7mm Hex
3.0	Metacarpal	3.9	2.0	3.0	25	70	1.0	2.0	T8 Torx
3.5	Metacarpal, Ulna	4.8	2.4	3.5	25	70	1.2	2.4	T10 Torx
4.0	Metacarpal, Ulna	5.5	2.75	4.0	35	80	1.2	2.75	T15 Torx
5.0	Ulna	6.9	3.5	5.0	35	90	1.6	3.5	T20 Torx

Note: Measurements are shown in mm

